



Measuring behavior, perception, attention and mental state

New sensors, tools and analytics support data-driven training and assessment

Lucas P.J.J. Noldus, Ph.D.

NextGen Training Technologies – Industrial Conference

13 June 2022

Company information

Noldus Information Technology

- Developer of software, hardware and integrated solutions for measurement and analysis of human behavior
- Founded in 1989
- Headquarters in Wageningen, The Netherlands
- Offices across Europe, North America, Asia
- 165 employees
- Clients: universities, research institutes, industry (automotive, aerospace, maritime, tech companies, etc.)



Human factors and interactive systems

Background and scope

- **Human factors:** application of **psychological** and **physiological** principles to the engineering and design of products, processes, and systems
- **Goal:** better understanding how humans can work in harmony with technology
- **Application:** system design, procedures, **training** to reduce human error, increase productivity, enhance safety, optimize comfort when operating an interactive system
- **Examples:** aircraft, car, train, ship, control room, or any other mission-critical system

Human factors approach to training and assessment

Processes to be addressed

- **Behavior**
- **Perception**
- **Attention**
- **Situational awareness**
- **Mental state**
- **Task performance**



Data-driven training and assessment

Evolution of assessment methods



Observation

Note taking
Qualitative analysis



Measurement

One or a few modalities
Single subject
Post-session analysis



Complex measurements

Multiple modalities
Multiple subjects
Real-time analysis and feedback

Human-system interaction

Measurement and analysis challenges

- What to measure?
- Which sensors to use?
- Synchronizing data streams
- Integrating data
- Analyzing data
- Obtaining insights



The Observer XT

Measuring human-system interaction

Technical trends

- Sensors become smaller, less expensive
- Sensor resolution goes up (e.g. video)
- Sensors become less invasive
- More computing power, cheaper storage
- Data integration: from offline to real-time
- Deep learning: automated pattern recognition

Measuring behavior, perception, attention and mental state

Eye tracking



Tobii Pro Spectrum



Tobii Pro Nano

Screen-based eye tracker

Contact-free
Fixed to one screen



Tobii Pro Glasses 3

Eye tracking glasses

Use anywhere
Operator free to move



Smart Eye Pro

Multi-camera eye tracker

Largest tracking distance,
motion box, view angle

Measuring behavior, perception, attention and mental state

Physiological sensing

- Cardiac activity - ECG (HR, HRV) or Pulse (HR)
- Skin conductance - GSR
- Muscle tension - EMG
- Respiration rate
- Brain activity – EEG, fNIRS



GSR – Galvanic skin response



Pulse rate



EMG



Multi-sensor wristbands



fNIRS



EEG



EEG

Measuring behavior, perception, attention and mental state

Facial imaging

Basic expressions

- Happy
- Sad
- Angry
- Surprised
- Disgusted
- Scared
- Contempt
- "Neutral"

Affective states

- Interest
- Confusion
- Boredom

Heart rate


- Remote PPG



FaceReader

Facial imaging

Facial action units




- FaceReader recognizes 20 facial action units
- See them in action on <https://www.noldus.com/facereader/facial-action-units>

Facial imaging

Beyond basic expressions


- Interest
- Boredom
- Confusion
- Attention
- Blink (AU45)
- Smiling
- Laughing
- Leaning Backward
- Leaning Forward
- Head turn left (AU51)
- Head turn right (AU52)
- Head up (AU53)
- Head down (AU54)



Highly relevant for human-system interaction

Facial imaging

Head orientation and gaze direction



Head orientation: recorded in three angles (X, Y and Z direction).


Detect nodding (drowsiness), turning (distraction)

Gaze direction: left-up, left, left-down, up, forward, down, right-up, right, right-down. Average error of ~7 degrees (~5 degrees or less for frontal faces~).

Detect areas of visual focus (not as accurate as an eye tracker)

Remote photoplethysmography

Reading your heart rate from a distance



- Blood pressure pulses → changes in blood volume → small changes in skin color, in sync with heart rate
- Software computes heart rate and heart rate variability
- Measure stress level, cognitive workload



Facial imaging

Test in real-world driving scenario




Daylight camera setup



Infrared camera setup

Facial imaging

Usage in car



- Driving context: Day
- Video camera: Daylight camera
- What this video tells us: FaceReader can detect driver drowsiness (eyes closed)

Facial imaging Usage in car



- Driving context: Day
- Video camera: Daylight camera
- What this video tells us: FaceReader can detect driver drowsiness (forward head rotation)

Facial imaging


Usage in car



- Driving context: Day
- Video camera: IR camera
- What this video tells us: FaceReader can detect the facial expressions “happy” and “surprised” in IR images captured in daylight

Facial imaging


Usage in car



- Driving context: Day
- Video camera: IR camera
- What this video tells us: FaceReader can detect head rotation (looking in external mirror) in IR images

Training operators, pilots and drivers

A multidimensional challenge




Requirements for assessment system

- Flexible: support many different scenarios, static/dynamic context, high/low workload
- Multimodal: different sensors
- Multiple subjects: interaction
- Real-time analysis and feedback
- Ease of use for the trainer

How to approach this challenge?

Context of the behavior

What is happening around the subject? Which factors are affecting him/her?



Measures representing behavior

Which data can be recorded that helps measure and understand the operator's actions?

Insights (also during the session)

How can we minimize the delay in a live feedback loop? Or how can we reduce time interpreting data?

Multimodal measurement

Concept validation

- Industry: Maritime (offshore)
- Use case: Operator in high-risk environment
- Simulation scenario: Tug boat captain has to go to an oil rig and position his barge to allow people to be lowered from the rig onto the deck of the barge.
- Measurement and analysis goals:
 - Where is he looking? Behavior (wearable eye tracking)
 - What is he looking at? Context (automated object detection)
 - Can he manage the task? Insight (cognitive load analysis)
- Technical challenge: automatic object / AOI labeling



The concept

Real-time contextual workload assessment by combining **attention**, using automated object labeling and gaze mapping, and **mental effort**, based on pupil-size change.

- Single sensor (head-mounted Eye Tracking)
- Automatic compensation of light intensity changes
- Real-time data processing

The Observer XT 14 - Marin2

File Edit View Setup Observe Analyze Help

Visualization - video1_detections_20

Project Explorer

Start episode Create episode selection

View Settings

Playback Control

02:26:22

1/2 x

Video

video1_detections_20.mp4

Automatically detect and identify displays and controls


Relative Time

02:26:22 (minutes:ff) 30 01:30:00 01:35:00 01:40:00 01:45:00 01:50:00 01:55:00 02:00:00 02:05:00 02:10:00 02:15:00 02:20:00 02:25:00 02:30:00 02:35:00 02:40:00 02:45:00 02:50:00 02:55:00 03:00:00 03:05:00 03:10:00 03:15:00 03:20:00 03:25:00

Results

Marin Objects

- outside
- telephone
- thruster
- rudder
- monitor:right_is...
- monitor:ship
- monitor:left_rader
- gyroscopes
- monitor:other
- overhead_sign
- monitor:dual
- monitor:map



The Observer XT 14 - Marin2

File Edit View Setup Observe Analyze Help

Visualization - video1_detections_20

Start episode Create episode selection View Settings

Playback Control

02:26:24 1/2 x

Video

video1_detections_20.mp4

Keep track of objects while operator moves around


Relative Time

02:26:24 (m:ss,ff) 00:01:30:00 01:35:00 01:40:00 01:45:00 01:50:00 01:55:00 02:00:00 02:05:00 02:10:00 02:15:00 02:20:00 02:25:00 02:30:00 02:35:00 02:40:00 02:45:00 02:50:00 02:55:00 03:00:00 03:05:00 03:10:00 03:15:00 03:20:00 03:25:00

Marin Objects

- outside
- telephone
- thruster
- rudder
- monitor:right_ra...
- monitor:ship
- monitor:left_radar
- gyroscopes
- monitor:other
- overhead_sign
- monitor:dual
- monitor:map

Results

The image shows a screenshot of 'The Observer XT 14' software. The main window displays a video titled 'video1_detections_20.mp4' showing a ship's interior with various objects highlighted by colored bounding boxes (green, blue, purple, orange). A white text box is overlaid on the left side of the video, containing the text 'Keep track of objects while operator moves around'. Below the video is a 'Results' panel with a 'Relative Time' axis from 02:26:24 to 03:25:00. This panel contains a legend for 'Marin Objects' and a corresponding activity timeline. The legend lists objects such as 'outside', 'telephone', 'thruster', 'rudder', 'monitor:right_ra...', 'monitor:ship', 'monitor:left_radar', 'gyroscopes', 'monitor:other', 'overhead_sign', 'monitor:dual', and 'monitor:map'. The timeline shows colored bars representing the duration of each object's activity.

The Observer XT 14 - Marin2

File Edit View Setup Observe Analyze Help

Visualization - videoT_detections_20

Project Explorer

Start episode Create episode selection View Settings

Playback Control

02:26:22 1/2 x

Video

videoT_detections_20.mp4

Automatically link gaze to objects

Relative Time

02:26:22 (minutes:ff) 30 01:30:00 01:35:00 01:40:00 01:45:00 01:50:00 01:55:00 02:00:00 02:05:00 02:10:00 02:15:00 02:20:00 02:25:00 02:30:00 02:35:00 02:40:00 02:45:00 02:50:00 02:55:00 03:00:00 03:05:00 03:10:00 03:15:00 03:20:00 03:25:00

Marin Objects

- outside
- telephone
- thruster
- rudder
- monitor:right_is...
- monitor:ship
- monitor:left_rader
- gyroscopes
- monitor:other
- overhead_sign
- monitor:dual
- monitor:map

Results

The Observer XT 14 - Marin2

File Edit View Setup Observe Analyze Help

Visualization - video1_detections_20

Project Explorer

Start episode Create episode selection

View Settings

Playback Control

02:26:22

1/2 x

Video

video1_detections_20.mp4

Generate a time line of viewed objects


Relative Time

02:26:22 (minutes:ff) 30 01:30:00 01:35:00 01:40:00 01:45:00 01:50:00 01:55:00 02:00:00 02:05:00 02:10:00 02:15:00 02:20:00 02:25:00 02:30:00 02:35:00 02:40:00 02:45:00 02:50:00 02:55:00 03:00:00 03:05:00 03:10:00 03:15:00 03:20:00 03:25:00

Marin Objects

- outside
- telephone
- thruster
- rudder
- monitor:right_is...
- monitor:ship
- monitor:left_rader
- gyroscopes
- monitor:other
- overhead_sign
- monitor:dual
- monitor:map


Results


The image shows a screenshot of 'The Observer XT 14' software. At the top, there's a menu bar with 'File', 'Edit', 'View', 'Setup', 'Observe', 'Analyze', and 'Help'. Below that is a toolbar with icons for 'Start episode' and 'Create episode selection'. A 'Playback Control' bar is visible with a progress slider and playback buttons. The main area is a video player showing a dark scene with several objects highlighted by colored bounding boxes (green, blue, yellow, orange). A white text box with a purple arrow points to the video player, containing the text 'Generate a time line of viewed objects'. Below the video player is a 'Results' section with a 'Relative Time' axis ranging from 02:26:22 to 03:25:00. This axis is populated with a timeline of colored bars corresponding to the objects in the legend. The legend lists 'Marin Objects' with color-coded squares: outside (green), telephone (blue), thruster (purple), rudder (cyan), monitor:right_is... (magenta), monitor:ship (red), monitor:left_rader (orange), gyroscopes (yellow), monitor:other (light green), overhead_sign (dark red), monitor:dual (light blue), and monitor:map (dark green).

- Based on pupil dilation
- Real-time and fully automated
- Automatic compensation of light intensity changes

Live cognitive load assessment


Identify peaks in workload





1x


Videos



Multimodal measurement and analysis

Preliminary results: attention


- 3 different levels of complexity:
 - › Easy: low wind & current
 - › Medium: average wind & current
 - › Hard: strong wind & current
- 3-stage task:
 - › Approach
 - › Maneuvering: positing the deck of the barge below the basket
 - › Basket lowering: keep the deck at the right position & make sure the timing of the basket lowering is in sync with the waves



Multimodal measurement and analysis

Preliminary results: Cognitive load


- Cognitive load derived from pupil dilation
- Corrected for changes in light exhibited by the displays



Where do we go from here?

Context of the behavior

What is happening around the subject? Which factors are affecting him/her?



Measures representing behavior


Which data can be recorded that helps measure and understand the operator's actions?

Insights (also during the session)

How can we minimize the delay in a live feedback loop? Or how can we reduce time interpreting data?


Human factors research and training

Traditional measurement setup



New: networked research and training platform

Noldus Hub









Noldus Hub

History Settings

testhub.noldus.com/projects/4880e900f-21aa0-4f51-8b30-bef6832631abf/inventory/50511/area-9c56-4303-af7a-bc1f63c05e51/record

DevOps TSD | Intranet Replianc Verlat Synges HR Noldus Hub Customer Journey... Dazelle Agir Projecten

Demo

NOLDUS HUB

Projects ▶


Setup

- Devices & Locations
- Members


Record

- Location overview
 - Location 1
 - Location 2

Screen



Webcam



Timeline

Select data

Heart rate	59 BPM
Interbeat interval	1016 ms
Tonic skin conductance	μ S 0.77

Record Session


CALIBRATE

T_3_09042022.apj

Show all X

Platform for data-driven training and assessment

Wrapping up: key benefits



Noldus Hub: work in progress

How can we work together?

- 
- Step 3
Let's design the approach
 - Step 2
Share your needs
 - Step 1
What is your challenge?

Thank you



Lucas P.J.J. Noldus, Ph.D.

lucas.noldus@noldus.nl | [linkedin.com/in/lucasnoldus](https://www.linkedin.com/in/lucasnoldus) | www.noldus.com



Noldus Information Technology



Radboud University
Nijmegen